



# Certificate of Traceability

## Manufacturer Certification

This certifies that *Tosoh Corporation*, using *1.6 min Variant Analysis Mode Ver. 5.24F on G8* has participated in and successfully completed the NGSP certification for manufacturers and is traceable to the **Diabetes Control and Complications Trial** Reference method. The comparison was performed with: **University of Missouri SRL#12**

The system evaluated was:

Instrument: <b>G8</b>	Calibrator Lot: <b>ZS3002</b>	Column Lot: <b>T</b>
Reagent Lot: <b>T8-101K, T8-203K, T8-301K, HW-43K</b>	Calibrator Assigned Values: <b>5.87%, 10.59%</b>	

Date of Certification: June 1, 2024

Certification Expires: June 1, 2025

A handwritten signature in black ink, appearing to read 'David Searles', is written over a horizontal line.

NGSP Steering Committee Chair

A handwritten signature in black ink, appearing to read 'Randie R Little PhD', is written over a horizontal line.

NGSP Network Coordinator

A handwritten signature in black ink, appearing to read 'Sharon Amalf', is written over a horizontal line.

SRL director/ supervisor