



Certificate of Traceability

Manufacturer Certification

This certifies that **Tosoh Corporation**, using **1.6 min Variant Analysis Mode Ver.5.24 on HLC-723G8** has participated in and successfully completed the NGSP certification for manufacturers and is traceable to the **Diabetes Control and Complications Trial** Reference method. The comparison was performed with: **University of Missouri SRL#9**

The system evaluated was:

Instrument: HLC-723G8	Calibrator Lot: ZS0002	Column Lot #: S
Reagent Lot: S8-101F, S8-201F, S8-301F, HW-19F	Calibrator Assigned Values: 5.95%, 10.31%	Ver.5.24

Date of Certification: June 1, 2021

Certification Expires: June 1, 2022

NGSP Steering Committee Chair

NGSP Network Coordinator

SRL director/ supervisor