


501 RP+



Paperless
chromatogram
graph

Assists
with auto-
verification

User-defined
Hb reference
library

QC
Levey-
Jennings

Reagent
management

Comprehensive Laboratory Management for Hemoglobin A1c Testing

TOSOH BIOSCIENCE

www.tosohbioscience.us

For In Vitro Diagnostic Use; Rx Only



What is 501RP+?

A data management application software to aid in workflow efficiency for your laboratory from processing to reporting.

AUTO-VERIFICATION

- Complex results flagging
- Customizable result rules

CHROMATOGRAM

- Customizable result reports
- Detailed peaks on graphs
- Overlay or comparison function

QUALITY CONTROL

- Centralized QC analysis
- Use of Levey-Jennings
- Employment of Westgard rules

GHb LIBRARY

- User-defined library
- Archive of abnormal hemoglobin

MAINTENANCE

- Embedded daily maintenance checklist
- Complete maintenance management
- On-board service records

TRACEABILITY

- Analysis trending of SA1c results
- Barcoded reagents and consumables
- On-board stability tracking and expiry

Available in Five Configurations of the Tosoh Automated Hemoglobin Analyzer Series



G11-90L*
G11-290L*



G8-90L
G8-290L



GX*

** Not FDA Cleared. Available outside US only*



"Tosoh Bioscience" is a registered trademark of Tosoh Corporation

Tosoh Bioscience, Inc.
6000 Shoreline Court, Suite 101
South San Francisco, CA 94080
Tel: (800) 248-6764
Fax: (650) 615-0415
www.tosohbioscience.us